Geistlich Bio-Gide® Compressed

ENGLISH

Geistlich Bio-Gide® Compressed

Resorbable Bilayer Collagen Membrane for Guided Tissue and Bone Regeneration

Geistlich Bio-Gide® Compressed is a biocompatible, sterile resorbable bilayer collagen membrane for guided tissue and bone regeneration.

COMPOSITION

Geistlich Bio-Gide® Compressed is a pure collagen membrane obtained by standardized, controlled manufacturing processes.

The collagen is extracted from veterinary certified pigs and is carefully purified to avoid antigenic reactions.

Geistlich Bio-Gide® Compressed sterilized in double blisters by γ -irradiation.

Geistlich Bio-Gide® Compressed has a bilayer structure: The porous surface – facing the bone – will allow the ingrowth of the bone forming cells. The dense surface – facing the soft tissue – will prevent the ingrowth of fibrous tissue into the bone defect. The membrane is made of collagen without further cross-linking.

PROPERTIES / ACTION

The low antigenicity and excellent biocompatibility support the use of Geistlich Bio-Gide® Compressed in dental surgery. Animal studies have shown that the collagen is resorbed in bone cavities within 24 weeks. The fibrous microstructure of Geistlich Bio-Gide® Compressed is hydrophilic. The membrane retains its structural integrity even when wet. Fixation by sutures or pins is possible. Geistlich Bio-Gide® Compressed has a smooth dense surface which facilitates cutting the membrane. The smooth, dense surface also prolongs the water uptake, application is facilitated. The coherent collagen fibres swell and form a unified basic tissue structure. As a result, the adaption to the bone wall and the complete closure of the bony defects is easily achieved. Inflammatory reactions have not been observed, but due to the special physical properties and elongated resorption time a certain inflammatory reaction after placing the membrane cannot be excluded.

INDICATIONS FOR USE

Geistlich Bio-Gide® Compressed is indicated for:

- augmentation around implants placed in immediate extraction sockets.
- augmentation around implants placed in delayed extraction sockets.
- localized ridge augmentation for later implantation.
- alveolar ridge reconstruction for prosthetic treatment.
- filling of bone defects after root resection, cystectomy, removal of retained teeth.
- guided bone regeneration in dehiscence defects.
- guided tissue regeneration procedures in periodontal defects.

CONTRAINDICATIONS

Geistlich Bio-Gide® Compressed should not be placed where active infection exists. Before placement, the surgeon should be confident that any active or recent infection has been properly treated.

Geistlich Bio-Gide® Compressed should not be used in patients with a known allergy to collagen.

DIRECTIONS FOR USE

The general principles of sterile handling and patient medication must be followed when using Geistlich Bio-Gide® Compressed:

- The bone defect is exposed by a mucoperiosteal flap and basic surgical procedures (e.g., curettage) are undertaken.
- The bone defect is filled with a space-making material, e.g., autogenous bone, Geistlich Bio-Oss®, Geistlich Bio-Oss Collagen®, etc. The defect should not be overfilled.
- The size of Geistlich Bio-Gide® Compressed is adapted by scissors according to the exposed defect. As a template a sterile aluminum foil may be used. The membrane should overlap the walls of the defect by at least two millimeters to allow complete bone contact and to prevent gingival connective tissue invasion below the material.
- The smooth surface (marked at the membrane margin with «UP») should face the soft tissue, and the rough side should face the bone.
- Geistlich Bio-Gide® Compressed is applied over the defect without further treatment and held in place with moderate pressure.
- The period of time necessary to apply pressure will vary with the degree of bleeding.
- Adhesion to the bone surface is facilitated by the hydrophilic properties of the membrane.
- Complete penetration of the membrane by blood and exudate allows excellent adaptation and adhesion of the membrane to bony structures.
- Salivary and other contamination to the material and surgical site should be minimized to avoid bacterial contamination.
- Due to the high tensile strength of this membrane, fixation is possible. Fixation may be indicated to avoid membrane displacement due to loading or mobilization.
- The mucoperiosteal flap is sutured over the collagen membrane (e.g., single sutures and deep mattress sutures).
- Complete wound closure should be attempted to avoid accelerated resorption due to membrane exposure, though clinical experience demonstrates satisfactory healing in most cases when complete wound closure is not possible.
- In order to avoid the formation of a long junctional epithelium when treating periodontal defects, the coronal portion of Geistlich Bio-Gide® Compressed should be adapted closely to the tooth (e.g. with additional fixation using suture material).

POST-OPERATIVE CARE

- The patient should be monitored closely.
- If the membrane becomes exposed, the dehiscence usually heals by itself within several weeks.
- Membrane removal is usually not necessary. However, to minimize bacterial contamination rinsing with bactericidal solutions is recommended.
- In the event that early membrane removal is necessary, the tissues adjacent to the membrane should be anesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual membrane. Following careful reflection of the surrounding tissue, the remaining portion of the membrane can be excised and the area curetted to remove any inflamed or infected tissue.
- To allow for undisturbed bone regeneration underneath the membrane surgical reentry should not be done before 4 to 6 months postoperatively.

PRECAUTIONS

Geistlich Bio-Gide® Compressed should only be used by trained dentists or oral surgeons. In case of membrane exposure during the healing phase the resorption time may be accelerated. A basic requirement for successful periodontal treatment includes eradicating the underlying bacterial infection as well as adequate oral hygiene. Therefore, prior to surgical intervention, patients must receive a hygiene phase of treatment, consisting of oral hygiene instructions, scaling and root planning, and occlusal adjustment when indicated. A postoperative maintenance phase can help to ensure long-term therapeutic success. If endosseous implants are involved, the membrane should be used only in combination with a stable implant and not in lieu of achieving primary implant stability. The long-term safety and effectiveness of the membrane in conjunction with bone-filling materials has not yet been established. Absolute stability of the membrane is important for guided bone regeneration and is a vital condition for therapeutic success, and the smallest movement on the tissue underneath is to be avoided. Effect on pregnant or lactating women is not known.

Geistlich Bio-Gide® Compressed should be used with special caution in patients with

- Acute or chronic infection (e.g., osteomyelitis) at the surgical site
- Uncontrolled metabolic diseases (e.g., diabetes, osteomalacia, thyroid disorder)
- Prolonged corticosteroid therapy
- Autoimmune diseases
- Radiotherapy
- Heavy smoking

ADVERSE REACTIONS

As Geistlich Bio-Gide® Compressed is a collagen product allergic reactions may not be totally excluded. Possible complications which may occur with any surgery include swelling at the surgical site, flap sloughing, bleeding, dehiscence, hematoma, increased sensitivity and pain, bone loss, redness, and local inflammation.

STORAGE AND HANDLING

If the sterile packaging is damaged or opened, the product must not be used. The contents of the double blister is designed for single use only.

Do not resterilize Geistlich Bio-Gide® Compressed. Resterilization is not possible. Discard any unused material. Do not use after expiry date.

Store at controlled room temperature 15°–25° C (59°–77°F) and in a dry place.

The membrane should be handled using sterile gloves or sterile instruments.

PRESENTATION

Geistlich Bio-Gide® Compressed is packed in sterile double blisters.

Re-order #	Size
500363	1 Geistlich Bio-Gide® Compressed membrane, 13 x 25 mm
500374	1 Geistlich Bio-Gide® Compressed membrane, 20 x 30 mm

Distributed by

Geistlich Pharma North America Inc.
Princeton, New Jersey 08540
1-(855) 799-5500
www.geistlich-na.com

Manufacturer

Geistlich Pharma AG
CH-6110 Wolhusen
Switzerland
Made in Switzerland

CAUTION: Federal law restricts this device to sale by or on the order of a dentist or physician.
Package insert: www.geistlich-na.com

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